

Human ACV-A (Activin A) QuickTest ELISA Kit

Catalog Number: HTK-YJL-0155

• Basic Information

Synonyms: Inhibin Beta A Chain ELISA Kit, Activin Beta-A Chain ELISA Kit, Inhba ELISA Kit

Species: Human

UniProt ID: P08476

Sample Type: Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 2 hours

Range: 15.625-1000 pg/ml (1 ng=1 mIU)

Sensitivity: 9.4 pg/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with Activin A, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.