

## Rat DRD3 (Dopamine Receptor D3) ELISA Kit

Catalog Number: HTK-YJL-0174

### • Basic Information

Synonyms: DRD3 ELISA Kit, Dopamine D3 Receptor ELISA Kit, D(3) Dopamine Receptor ELISA Kit

Species: Rat

UniProt ID: P19020

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 4 hours

Range: 0.313-20 ng/ml

Sensitivity: 0.188 ng/ml

Intended Use: For research use only.



### • Specificity & Cut Off

Specificity: Specifically binds with Drd3, no obvious cross reaction with other analogues.

### • Specification

48T, 96T

### • Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnostech.com](mailto:info@altadiagnostech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA