

Rat E2 (Estradiol) ELISA Kit

Catlog Number: HTK-YJL-0177

• Basic Information

Synonyms: Estradiol ELISA Kit, E2 ELISA Kit

Species: Rat

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Competitive ELISA, Coated with Antigen

Detection Wavelength: OD450

Reaction Duration: 2 hours

Range: 12.5-800 pg/ml

Sensitivity: 7.5 pg/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with E2, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.