

Rat PI (Proinsulin) ELISA Kit

Catalog Number: HTK-YJL-0181

• Basic Information

Synonyms: PI ELISA Kit, Proinsulin ELISA Kit

Species: Rat

UniProt ID: P01322

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 4 hours

Range: 0.313-20 ng/ml

Sensitivity: 0.188 ng/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with PI, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.