

## Human DBH (Dopamine Beta-hydroxylase) ELISA Kit

Catalog Number: HTK-YJL-0182

### • Basic Information

Synonyms: Dopamine Beta-hydroxylase ELISA Kit, Dopamine Beta-monoxygenase ELISA Kit, Soluble Dopamine Beta-hydroxylase ELISA Kit, DBH ELISA Kit

Species: Human

UniProt ID: P09172

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 4 hours

Range: 1.563-100 ng/ml

Sensitivity: 0.938 ng/ml

Intended Use: For research use only.



### • Specificity & Cut Off

Specificity: Specifically binds with DBH, no obvious cross reaction with other analogues.

### • Specification

48T, 96T

### • Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnostech.com](mailto:info@altadiagnostech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA