

Hamster ADP/Acrp30 (Adiponectin) ELISA Kit

Catalog Number: HTK-YJL-0203

• Basic Information

Synonyms: Adiponectin ELISA Kit, ADIPOQ ELISA Kit, ACDC ELISA Kit, ACRP30 ELISA Kit, APM1 ELISA Kit, GBP28 ELISA Kit, 30 KDa Adipocyte Complement-related Protein ELISA Kit, Adipocyte Complement-related 30 KDa Protein ELISA Kit, ACRP30 ELISA Kit, Adipocyte, C1q And Collagen Domain-containing Protein ELISA Kit, Adipose Most Abundant Gene Transcript 1 Protein ELISA Kit

Species: Hamster

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 4 hours

Range: 15.625-1000 pg/ml

Sensitivity: 9.375 pg/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with ADP/Acrp30, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA